

Smith+Nephew

A closer look at the use of placental membranes for wound closure

Charles B. Parks, DPM, FACFAS



Table of contents:

Introduction	3
Chronic wounds in the United States: Prevalence, risk factors, and best practices	3
Why use placental tissues in the treatment of chronic wounds?	4
Why use GRAFIX PL [®] Lyopreserved and GRAFIX [®] Cryopreserved Placental Membranes in the treatment of chronic wounds?	5
What is the evidence for using GRAFIX PL and GRAFIX Membranes?.....	6
Why use STRAVIX [®] Cryopreserved Umbilical Tissue in the treatment of chronic wounds?	7
What is the evidence for using STRAVIX [®] Tissue?	8
Case studies	9
Conclusion.....	16



Charles B. Parks, DPM, FACFAS
**Chief of Podiatric Surgery, Zuckerberg San Francisco
General Hospital and Trauma Center**
**Assistant Clinical Professor, University of California,
San Francisco**
**Assistant Clinical Professor, California School of
Podiatric Medicine**

Charles B. Parks, DPM, is one of the treating physicians at Zuckerberg San Francisco General Hospital and Trauma Center (SFGH)/University of California, San Francisco (UCSF) in the Department of Orthopedics and specializes in treating lower extremity wounds in patients at high risk of amputation, in patients requiring forefoot reconstruction, and in sports medicine. Dr. Parks completed his three-year surgical residency in podiatric surgery at the Veterans Hospital in San Francisco (VASF) where he worked at many different surgery centers and hospitals in the Bay Area, including VASF, UCSF, and SFGH. In 2016, he joined the UCSF faculty as an Assistant Clinical Professor where he has continued his research and teaching while also holding an Assistant Clinical Professor role at the California School of Podiatric Medicine. In 2018, he was named Chief of Podiatric Surgery at Zuckerberg San Francisco General Hospital. Dr. Parks is a member of numerous professional associations and is double board certified by the American College of Foot and Ankle Surgeons and the American Board of Podiatric Medicine.

Introduction

Chronic wounds, sometimes referred to as stalled wounds, are characterized by a failure to reduce in size by 40% in 30 days. In 2019, there were 8.2 million people in the United States with chronic wounds.¹ Annually, about 1.5 million Americans experience a diabetic foot ulcer and 1 million a venous leg ulcer.²

Such wounds can greatly benefit from the application of placental tissues, such as GRAFIX PL[◊] and GRAFIX[◊] Membranes and STRAVIX PL[◊] and STRAVIX[◊] Tissue by Smith+Nephew. GRAFIX PL, GRAFIX Membranes, STRAVIX PL, and STRAVIX Tissue all contain native living cells, growth factors, and an extracellular matrix. These products also promote re-epithelization and speed wound closure, even for the most complex of chronic wounds.

Chronic wounds in the United States: Prevalence, risk factors, and best practices

Lower extremity wounds tend to have poor outcomes: for patients with diabetic foot ulcers, there is a 40% recurrence rate within one year.³ This is particularly troubling because some 85% of all lower limb amputations are preceded by a diabetic foot ulceration.⁴⁻⁶ For patients with venous leg ulcers, there is a 70% recurrence rate within three months, and 66% of such wounds fail to heal with standard of care in 12 weeks.⁷⁻¹¹ Ultimately, each day a chronic wound remains open, the risk that the patient will develop an infection, be hospitalized, or undergo an amputation increases.

Chronic wounds negatively impact most patients who have them. In fact, 68% of people with chronic leg wounds reported that their life was negatively affected by the wound.¹² Negative effects of chronic wounds on patients include reduced quality of life, as well as reduced worker productivity that can lead to financial burden. Risk factors for developing chronic wounds include obesity, diabetes, smoking, hypertension, and immune deficiency.¹³⁻¹⁸

For chronic wounds to heal, the wound bed must be cared for according to the TIME principles. **T** stands for tissue; **I**, infection and/or inflammation; **M**, moisture imbalance; and **E**, edge of wound.¹⁹⁻²³

According to the TIME principles, wounds have the best chance at healing when:

- Non-viable tissue is debrided.
- Infection or inflammation is managed through antimicrobial agents or further debridement.
- Moisture balance is maintained.
- An intact epithelium is re-established, restoring skin function.

Re-establishing an intact epithelium (or re-epithelization) can be facilitated in chronic wounds by using placental tissues as a wound cover, wrap, or barrier.

Why use placental tissues in the treatment of chronic wounds?

Fresh placental tissues contain native living cells, such as epithelial cells, fibroblasts and mesenchymal stem cells, a collagen-rich extracellular matrix, growth factors, and other mediators.

Since the early twentieth century, placental tissues have been used as a biological dressing to care for wounds, treat burns, and reconstruct ocular surfaces. Although the thinner placental membranes (which consist of amnion and chorion) are utilized most often, there is an increasing interest in using umbilical cord tissue for surgical applications, particularly limb salvage cases, because it is thicker and more durable than placental membranes.


Donors of placental tissues are screened for most diseases.

Preserving placental tissues


Before the placental tissue can be used in a clinical setting, it must be preserved to make a shelf-stable commercial product. Most placental tissue products undergo traditional tissue lyophilization, which can kill native living cells, remove growth factors, and alter the structure of the extracellular matrix.²⁴

Fortunately, there have been advances in preservation techniques in recent years. Prestige Lyotechnology[®], for example, preserves the extracellular matrix, growth factors, and native living cells. This technology produces a product that can be stored at room temperature, thus eliminating the need for ultralow-temperature storage.

Prestige Lyotechnology, a type of lyopreservation, is the preservation method for GRAFIX PL[®] Membrane and STRAVIX PL[®] Tissue. Cryopreservation is the preservation method for GRAFIX[®] Membrane and STRAVIX[®] Tissue, both of which must be stored in a freezer until use. Both the cryopreservation and lyopreservation processes retain the placental tissue extracellular matrix, growth factors and other mediators, and native living cells.



**GRAFIX PL[®] Lyopreserved
Placental Membrane**



**STRAVIX[®] Cryopreserved
Umbilical Tissue**

Why use GRAFIX PL[◊] and GRAFIX[◊] Membranes to treat chronic wounds?


GRAFIX PL and GRAFIX Membranes are human placental membranes composed of an extracellular matrix, growth factors, and native living cells. They were developed as wound covers for a wide variety of acute and chronic wounds in head-to-toe locations. They naturally conform to complex anatomic locations and may be used over exposed bone, tendon, joint capsule, muscle, and hardware. Each lot is tested to ensure that native living cells are present across at least 70% of the membrane.

GRAFIX PL PRIME[◊] and GRAFIX PRIME[◊] Membranes are derived from amnion. They contain an epithelial layer with epithelial cells, as well as a stromal layer with fibroblasts and mesenchymal stem cells. GRAFIX CORE[◊] Membrane is derived from chorion. It contains a stromal layer with fibroblasts and mesenchymal stem cells, is trophoblast free, and contains no

epithelial cells. The three products differ only slightly in terms of their tensile strength and handling properties, but they are structurally and functionally equivalent in terms of efficacy and outcomes.

Both GRAFIX PL PRIME and GRAFIX PRIME Membranes come in different sizes, ranging from 1.5 cm x 2 cm to 5 cm x 5 cm. Both are also offered in 16 mm discs. GRAFIX CORE Membrane is available in 3 cm x 4 cm and 5 cm x 5 cm sizes. A variety of sizes allows for minimal waste and decreased costs to the health care system. GRAFIX PL PRIME Membrane is packaged between two mesh sheets that can be easily peeled apart to release the membrane, whereas GRAFIX PRIME and GRAFIX CORE Membranes are supplied on a plastic applicator.

Cross-section of placental membranes



What is the evidence for using GRAFIX PL^o and GRAFIX^o Membranes?

Robust evidence supports the effectiveness of GRAFIX PL and GRAFIX Membranes. In 2019, Ananian et al. published a retrospective, open-label, multicenter study of GRAFIX PL Membrane in the management of chronic wounds that had not responded to standard of care.²⁵ There were 78 patients in the study, with 98 wounds, and most of the patients had more than two comorbidities. The wounds included venous leg ulcers, diabetic foot ulcers, and surgical wounds, among other wound types, with an average wound area of 13.3 cm². The mean wound duration before treatment with GRAFIX PL Membrane was just under nine months. The authors showed that the application of GRAFIX PL Membrane caused 59% of non-healing wounds to achieve complete closure. It took median 63 days for wounds treated with GRAFIX PL Membrane to achieve closure, with median six graft applications.


104 such wounds.²⁷ The wounds were 12.7 cm² in size, with a wound duration before treatment of 3.9 months. According to these authors, 83.7% of the wounds—burns, pressure injuries, surgical and traumatic wounds, and necrotizing fasciitis, among others—completed closure within median 41 days, with a median of three grafts per wound.

Patient Demographics and Baseline Wound Characteristics


Patients	78
Wounds	98
Age (mean)	62.7 years
Patients with ≥ 2 comorbidities	65%
Wound area (mean)	12.7 cm²
Wound duration prior to treatment (mean)	8.9 months
Wound >12 months duration	20.4%

Most of the clinical studies conducted to analyze the impact of GRAFIX PL and GRAFIX Membranes on chronic wounds have tended to focus on diabetic foot and venous leg ulcers. For that reason, the recent publication of a study in *Wounds* is particularly welcome: it shows that non-diabetic and non-venous chronic wounds of the upper and lower extremities also respond well to GRAFIX Membranes.²⁷ Johnson et al. performed a study involving 92 patients with

Proportion of Patients Who Achieved Closure²⁵




Median Days to Closure²⁵⁻²⁶




There are many large clinical trials and studies that analyze the efficacy of GRAFIX PL and GRAFIX Membranes.^{25, 27-33} Altogether, according to these studies, the closure rates for wounds treated with the membranes range from 53% for the most complicated refractory venous leg ulcers to up to 83% for non-diabetic and non-venous chronic wounds.

Wound Closure by Location²⁷



Wound Closure by Wound Type²⁷



Why use STRAVIX^o Tissue in the treatment of chronic wounds?

Like GRAFIX PL and GRAFIX Membranes, STRAVIX Tissue is preserved to maintain the extracellular matrix, growth factors, and native living cells. It is conforming, adapting to injured tissue. STRAVIX Tissue is well suited for surgical soft tissue repair in the operating room such as limb salvage, and it is commonly used across many surgical specialties, with usage including wraps over tendon repair sites.

STRAVIX Tissue is composed of the umbilical amnion and Wharton’s jelly of human umbilical cord (which contains hyaluronic acid). To prevent an immune response in the patient receiving STRAVIX Tissue, the blood vessels in the umbilical cord are removed. Because the placement of the blood vessels within the umbilical cord varies from donor to donor, the appearance and thickness of the final product may vary. It is available in three sizes (2 cm x 2 cm, 2 cm x 4 cm, and 3 cm x 6 cm) and can be sewn or stapled together to make larger pieces or trimmed to make it smaller.

Patient Demographics and Baseline Wound Characteristics

Patients	10
Age	59.8 years
BMI (mean)	28.6 kg/m ²
NPWT Used	9/10 patients
Wound size	45.9 cm²

What’s the difference between GRAFIX PL and GRAFIX Membranes and STRAVIX Tissue? The primary difference is that STRAVIX Tissue is more durable. It is also thicker, ranging from 1 mm to 3 mm in thickness, and has greater tensile strength than GRAFIX PL and GRAFIX Membranes.

What is the evidence for using STRAVIX^o Tissue?

In 2018, Dr. McGinness investigated the use of STRAVIX Tissue to manage large, open soft tissue defects with gas gangrene.³⁴ The 10 patients chosen for the study had recently undergone surgery that resulted in a large, open defect with exposed structures, with wounds averaging 45.9 cm² in size.

All patients in the study had a history of diabetes, and many had coronary heart disease, congestive heart failure, hypertension, and other comorbidities. Because the patients had inadequate surrounding soft tissue or skin, surgical wound closure was not possible. The patients were not candidates for autologous skin grafts or flaps.

Before treatment with STRAVIX Tissue, the patients received aggressive debridement to remove infected, necrotic, and devascularized tissue. After irrigation by motorized pulse lavage, STRAVIX Tissue was fenestrated or sutured into the wound bed. The patients were administered intravenous antibiotics during their hospital stay, and nine of 10 of these patients received negative pressure wound therapy (NPWT) after surgery. It is important to note that there was no debridement performed after surgery.

Each patient achieved complete wound closure with one application of STRAVIX Tissue in 13.4 weeks on average. Remarkably, the mean length of these patients' hospitalization was only nine days. This is particularly significant compared to a large database review of over 11,000 patients with gas gangrene of the lower extremity that found the average length of hospital stay to be 31.1 days.³⁵

In 2020, Mullins et al. investigated the use of STRAVIX Tissue in the management of complex wounds with exposed bone, tendon, hardware, or other underlying structures.³⁶ The study was designed to be a quick pilot study to gather preliminary prospective data on a single application of STRAVIX Tissue. Ten patients with 12 wounds were treated over a four-week period. Wounds included venous leg ulcers, diabetic foot ulcers, pressure injuries, traumatic wounds, and one surgical wound. The average wound size was 16.5 cm², and the mean wound duration before treatment was 10 months. Eighty percent of the wounds achieved full granulation in four weeks or less. The time to complete granulation was about 17.5 days. In addition, 30% of the wounds achieved complete closure in four weeks or less.

These data show good progression to wound closure in a short period of time. If the patients had been followed for longer, it is anticipated that more patients would have achieved complete wound closure, which occurred in the study by McGinness et al. discussed earlier.³⁴ Moreover, the data in this study support the use of STRAVIX Tissue as a cover until the wound is fully granulated, after which the wound could be ready for a skin graft or flap at the physician's preference.



Case studies

Case study no. 1: Treatment of a chronic wound with GRAFIX PL^o Membrane

A 48-year-old woman presented with a 9 cm² chronic dorsal foot wound with a mixed fibrogranular wound base. The patient had a history of foot infection and previously underwent a partial first-ray resection. She had uncontrolled type 2 diabetes with glycosylated hemoglobin (HbA1C) of 13.5. She was also obese and had hypertension.

The wound was first irrigated, and all non-viable soft tissue was debrided. Then, after thorough excisional debridement, a GRAFIX PL Membrane was applied over the wound. After two applications of GRAFIX PL Membrane with serial debridement, the wound base was granular, and the wound was smaller. In total, there were six weekly applications of the GRAFIX PL Membrane. At six weeks, there was complete wound closure.



1. Initial Presentation



2. Application #2



3. Application #3



4. Application #4



5. Application #6



6. One week later



8. Closure achieved

Case study no. 2: Treatment of an amputation site with GRAFIX[®] Membrane and STRAVIX[®] Tissue

A 59-year-old man presented with forefoot gangrene with a surgical defect. He had type 2 diabetes, peripheral vascular disease, and end-stage renal disease requiring daily hemodialysis. Before amputation of the forefoot, he underwent an angioplasty to revascularize the leg.

After amputation, there was a 32 cm² dorsal foot wound, which was covered with STRAVIX Tissue. NPWT was applied for one week. At four weeks, the graft was still incorporating, as the wound continued to shrink in size.

At six weeks, the wound was fully granular with a little bit of the plantar fat pad exposed. The initial thought was to take the patient back to the operating room to debride the fat pad and perhaps regrant the wound. However, by the time of the surgery, most of the fat pad was granulated over. For that reason, minimal debridement was performed, and the fat pad was tacked in place. Then, GRAFIX PRIME[®] Membrane was placed over the wound.

Ten weeks after the STRAVIX Tissue application and one week after the GRAFIX PRIME Membrane application, the grafts were fully incorporated. The wound continued to decrease in size, and, at week 18, there was complete wound closure. In total, there was one STRAVIX Tissue application and one GRAFIX PRIME Membrane application.



1. Amputation site



2. Initial Amputation



3. Application of STRAVIX[®] Tissue



4. 12 days after STRAVIX Tissue application



5. 4 weeks after STRAVIX[®] Tissue placement



6. 6 week follow up



7. OR debridement with GRAFIX PRIME[®] membrane placement at 9 weeks



8. 10 weeks post GRAFIX PRIME Membrane



9. 12 weeks post GRAFIX PRIME Membrane



10. 14 weeks post GRAFIX PRIME Membrane



11. Closure at 18 weeks

Case study no. 3: Treatment of a postoperative infected wound with STRAVIX[®] Tissue

A 48-year-old man presented with a Lisfranc fracture that was treated with open reduction and internal fixation surgery. The patient was HIV positive and had progressive macular hypomelanosis. He also had bipolar disorder and polysubstance dependence.

After surgery, the patient did not attend post-surgery follow-up appointments. When the patient eventually returned for treatment, he presented with a deep-space infection that went down to the hardware and joints. He underwent irrigation and debridement to remove the non-viable tissue and hardware. He also received six weeks of intravenous antibiotics (nafcillin) to resolve the infection.

He ultimately had two wounds on his foot—one dorsal and the other medial—that went down to the joint and joint capsule. One application of STRAVIX Tissue was applied with a non-adherent dressing, and NPWT was applied for one week. At one month, the wounds were smaller. More importantly, the wounds were completely granular and filled in. At four months, both wounds were fully closed.



1. Preoperative



2. Procedure



**Initial presentation
STRAVIX® Tissue application**

1 week after application

**One-month status post-STRAVIX®
Tissue application**



**Two months status post-STRAVIX® Tissue
application**

**Three months status
post-STRAVIX® Tissue application**

**Four months -
Closure of both wounds**

Case study no. 4: Treatment of a posterior ankle wound with an exposed Achilles tendon with STRAVIX^o Tissue

A 59-year-old woman presented with a full-thickness posterior ankle ulcer with exposed Achilles tendon. The patient had type 2 diabetes, chronic limb-threatening ischemia, calciphylaxis, hypertension, and end-stage renal disease that required hemodialysis three times a week. The ankle ulcer was caused by pressure and calciphylaxis, and the soft tissue was affected with methicillin-resistant *Staphylococcus aureus* (MRSA).

After two sharp debridements were performed, the wound left the Achilles tendon exposed. The wound was covered in two pieces of STRAVIX Tissue that were stapled together, over which a non-adherent dressing was placed. STRAVIX Tissue was chosen because of the size and depth of the ulcer and because the patient had complex comorbidities known to affect wound healing negatively. After STRAVIX Tissue was applied, the patient underwent NPWT for one week. In addition, the patient underwent an angiographic procedure with interventional popliteal angioplasty for near-total occlusion of her left tibial-peroneal trunk. She also received two weeks of post-hemodialysis antibiotics (vancomycin).

After six weeks, the wound was smaller and granular, with some creeping granulation over the Achilles tendon. The wound continued to grow smaller, and, at 15 weeks, complete wound closure occurred. One year later, the wound was still closed, and she had good range of motion of the ankle joint.



1. Initial presentation



2. Follow up: s/p 2nd sharp debridement with 5 days of NPWT



3. Pre-application



4. Application of STRAVIX® Tissue
8 cm x 4 cm x 0.7 cm



5. 6 weeks after application



6. 7.5 weeks
6.2 cm x 2.1 cm x 0.2 cm



7. 8.5 weeks
5.5 cm x 2.1 cm x 0.2 cm



8. 9.5 weeks
Proximal: 1.5 cm x 1.5 cm x 0.1 cm
Distal: 2 cm x 2.2 cm x 0.1 cm



9. Closed at 15 weeks



10. Remains closed at 1 year

Case study no. 5: Treatment of an extensor hallucis longus rupture repair with STRAVIX® Tissue

A 42-year-old man presented with degeneration of the talonavicular joint after talonavicular fusion surgery with a dorsal spanning plate. Unfortunately, the plate ruptured his extensor hallucis longus tendon. By the time he was taken to the operating room, which was only one week after the rupture, there was a 4 cm gap from end to end.

After the extensor hallucis longus tendon repair, the repair and the tendon itself were wrapped in STRAVIX Tissue and sutured in place. The patient was kept non-weight bearing in a short leg cast for six weeks. At six weeks, he was transferred to a controlled ankle motion (CAM) boot. Full closure of the wound occurred eight weeks after surgery.



1. EHL Rupture Repair



2. Intra-op



5. Eight weeks post-op

Conclusion

A growing body of research shows that chronic wounds that are challenging to treat may respond well to placental tissues. There is a variety of skin substitute products available on the market, but, unfortunately, the quality of most of these products is greatly reduced during the preservation process.

Thanks to the advanced preservation methods pioneered by Smith+Nephew, GRAFIX PL[®] and GRAFIX[®] Membranes and STRAVIX[®] Tissue retain all components of fresh placental tissue, including native living cells, as well as growth factors and an intact extracellular matrix. Studies have proven that the Smith+Nephew suite of placental tissues hastens wound closure, thus improving patients' health and quality of life.

Studies have demonstrated that placental tissues such as those used in GRAFIX PL, GRAFIX Membranes, STRAVIX PL, and STRAVIX Tissue hasten wound closure.^{27,31} Achieving wound closure can help patients get back to their normal daily routines, improving quality of life.

References 1. Sen CK. Human wounds and its burden: an updated compendium of estimates. *Adv Wound Care (New Rochelle)*. 2019; 8(2): 39-48. 2. SmartTRAK™. <https://app.smarttrak.com/markets/qs/2f491c456352f41ac0a80237b1942c19> 3. Armstrong DG, Boulton AJM, Bus SA. Diabetic foot ulcers and their recurrence. *N Engl J Med*. 2017; 376(24): 2367-2375. 4. Frykberg RG, Zgonis T, Armstrong DG, et al. Diabetic foot disorders. A clinical practice guideline (2006 revision). *J Foot Ankle Surg*. 2006; 45(5): S1-S66. 5. Singh N, Armstrong DG, Lipsky BA. Preventing foot ulcers in patients with diabetes. *JAMA*. 2005; 293(2): 217-228. 6. Apelqvist J, Larsson J. What is the most effective way to reduce incidence of amputation in the diabetic foot? *Diabetes Metab Res Rev*. 2000; 16(S1): S75-S83. 7. Greer N, Foman NA, MacDonald R, et al. Advanced wound care therapies for nonhealing diabetic, venous, and arterial ulcers: a systematic review. *Ann Intern Med*. 2013; 159(8): 532-542. 8. Mostow EN, Haraway GD, Dalsing M, et al. Effectiveness of an extracellular matrix graft (OASIS Wound Matrix) in the treatment of chronic leg ulcers: a randomized clinical trial. *J Vasc Surg*. 2005; 41(5): 837-843. 9. Harding K, Sumner M, Cardinal M. A prospective, multicentre, randomised controlled study of human fibroblast-derived dermal substitute (Dermagraft) in patients with venous leg ulcers. *Int Wound J*. 2013; 10(2): 132-137. 10. Bianchi C, Cazzell S, Vayser D, et al. A multicentre randomised controlled trial evaluating the efficacy of dehydrated human amnion/chorion membrane (EpiFix®) allograft for the treatment of venous leg ulcers. *Int Wound J*. 2018; 15(1): 114-122. 11. Franks PJ, Barker J, Collier M, et al. Management of patients with venous leg ulcers: challenges and current best practice. *J Wound Care*. 2016; 25 Suppl 6: S1-S67. 12. Phillips T, Stanton B, Provan A, et al. A study of the impact of leg ulcers on quality of life: financial, social, and psychologic implications. *J Am Acad Dermatol*. 1994; 31(1): 49-53. 13. Guo S, DiPietro LA. Factors affecting wound healing. *J Dent Res*. 2010; 89(3): 219-229. 14. Pocock ES, Alsaigh T, Mazor R, et al. Cellular and molecular basis of venous insufficiency. *Vasc Cell*. 2014; 6(1): 24. 15. Strbo N, Yin N, Stojadinovic O, et al. Innate and adaptive immune responses in wound epithelialization. *Adv Wound Care (New Rochelle)*. 2014; 3(7): 492-501. 16. Sheehan PJ, Jones P, Caselli A, et al. Percent change in wound area of diabetic foot ulcers over a 4-week period is a robust predictor of complete healing in a 12-week prospective trial. *Diabetes Care*. 2003; 26(6): 1879-1882. 17. van Rijswijk L, Polansky M. Predictors of time to healing deep pressure ulcers. *Ostomy Wound Manage*. 1994; 40(8): 40-48. 18. Margolis DJ, Allen-Taylor L, Hoffstad O, et al. The accuracy of venous leg ulcer prognostic models in a wound care system. *Wound Repair Regen*. 2004; 12(2): 163-168. 19. Schultz GS, Sibbald RG, Falanga V, et al. Wound bed preparation: a systematic approach to wound management. *Wound Repair Regen*. 2003; 11: S1-S8. 20. Leaper DJ, Schultz G, Carville K, et al. Extending the TIME concept: what have we learned in the past 10 years? *Int Wound J*. 2012; 9(Suppl 2): 1-19. 21. Edmonds M, Foster AVM, Vowden P. Wound bed preparation for diabetic foot ulcers. In: *Wound Bed Preparation in Practice*. London, United Kingdom: MEP Ltd.; 2004: 6-11. 22. Kirsner RS. The standard of care for evaluation and treatment of diabetic foot ulcers. <http://www.barry.edu/includes/docs/continuing-medical-education/diabetic.pdf>. Accessed December 11, 2020. 23. Adeshara KA, Diwan AG, Tupe RS. Diabetes and complications: cellular signaling pathways, current understanding and targeted therapies. *Curr Drug Targets*. 2016; 17(11): 1309-1328. 24. Dhall S, Sathyamoorthy M, Kuang JQ, et al. Properties of viable lyopreserved amnion are equivalent to viable cryopreserved amnion with the convenience of ambient storage. *PLOS One*. 2018; 13(10): e0204060. 25. Ananian CE, Davis RD, Johnson EL, et al. Wound closure outcomes suggest clinical equivalency between lyopreserved and cryopreserved placental membranes containing viable cells. *Adv Wound Care (New Rochelle)*. 2019; 8(11): 546-554. 26. Data on file at Osiris. 27. Johnson EL, Saunders MC, Thote T, et al. Cryopreserved placental membranes containing viable cells result in high closure rate of non-healing upper and lower extremity wounds of non-diabetic and non-venous pathophysiology. *Wounds*. Forthcoming. 28. Lavery LA, Fulmer J, Shebetka KA, et al. The efficacy and safety of Grafix® for the treatment of chronic diabetic foot ulcers: results of a multi-centre, controlled, randomised, blinded, clinical trial. *Int Wound J*. 2014; 11(5): 554-560. 29. Frykberg RG, Gibbons GW, Walters JL, et al. A prospective, multicentre, open-label, single-arm clinical trial for treatment of chronic complex diabetic foot wounds with exposed tendon and/or bone: positive clinical outcomes of viable cryopreserved human placental membrane. *Int Wound J*. 2016; doi: 10.1111/iwj.12649. 30. Raspovic KM, Wukich DK, Naiman DQ, Lavery LA, Kirsner RS, Kim PJ, Steinberg JS, Attinger CE, Danilkovitch A. Effectiveness of viable cryopreserved placental membranes for management of diabetic foot ulcers in a real-world setting. *Wound Repair Regen*. 2018 Apr 23. doi: 10.1111/wrr.12635. 31. Farivar BS, Toursavadkoshi S, Monahan TS, et al. Prospective study of cryopreserved placental tissue wound matrix in the management of chronic venous leg ulcers. *J Vasc Surg Venous Lymphat Disord*. 2019; 7(2): 228-233. 32. Regulski M, Jacobstein DA, Petranto RD, et al. A retrospective analysis of a human cellular repair matrix for the treatment of chronic wounds. *Ostomy Wound Manage*. 2013; 59(12): 38-43. 33. Johnson E, Marshall J, Michael G. A comparative outcomes analysis evaluating clinical effectiveness in two different human placental membrane products for wound management. *Wound Repair Regen*. 2016; doi: 10.1111/wrr.12503. 34. McGinness K, Kurtz Phelan DH. Use of viable cryopreserved umbilical tissue for soft tissue defects in patients with gas gangrene: a case series. *Wounds*. 2018; 30(4): 90-95. 35. Fincke BG, Miller DR, Turpin R. A classification of diabetic foot infections using ICD-9-CM codes: Application to a large computerized medical database. *BMC Health Serv Res*. 2010; 10: 192. 36. Mullins RF, Hassan Z, Homsombath B, et al. Prospective, single-center, open-label, pilot study using cryopreserved umbilical tissue containing viable cells in the treatment of complex acute and chronic wounds. *Wounds*. 2020; 32(8): 221-227.

Advanced Wound Management

Smith+Nephew, Inc.
Fort Worth, TX 76109 USA

www.smith-nephew.com

www.grafixpl.com

GRAFIX and STRAVIX Customer Care Center

T 888-674-9551
F 443-283-4419

◇Trademark of Smith+Nephew
Osiris, GRAFIX, GRAFIX PL, GRAFIX
PRIME, GRAFIX CORE, GRAFIX PL
PRIME, and STRAVIX are trademarks
of Osiris Therapeutics, Inc., a wholly
owned subsidiary of Smith+Nephew
Consolidated, Inc.

© 2021 Smith+Nephew

GRCE51-22545-0321